

## カテーテルサイズ換算表

EG:English Gauge イングリッシュサイズ

Fr:French Gauge フレンチサイズ

EG	外径(mm)	Fr
-	0.3	1
00	0.5	-
-	0.7	2
0	1.0	3
-	1.3	4
1	1.5	-
-	1.7	5
2	2.0	6
-	2.3	7
3	2.5	-
-	2.7	8
4	3.0	9
-	3.3	10
5	3.5	-
-	3.7	11
6	4.0	12
-	4.3	13
7	4.5	-
-	4.7	14
8	5.0	15
-	5.3	16
9	5.5	-
-	5.7	17
10	6.0	18
-	6.3	19
11	6.5	-
-	6.7	20

EG	外径(mm)	Fr
12	7.0	21
-	7.3	22
13	7.5	-
-	7.7	23
14	8.0	24
-	8.3	25
15	8.5	-
-	8.7	26
16	9.0	27
-	9.3	28
17	9.5	-
-	9.7	29
18	10.0	30
-	10.3	31
19	10.5	-
-	10.7	32
20	11.0	33
-	11.3	34
21	11.5	-
-	11.7	35
22	12.0	36
-	12.3	37
23	12.5	-
-	12.7	38
24	13.0	39
-	13.3	40
25	13.5	-

EG	外径(mm)	Fr
-	13.7	41
26	14.0	42
-	14.3	43
27	14.5	-
-	14.7	44
28	15.0	45
-	15.3	46
29	15.5	-
-	15.7	47
30	16.0	48
-	16.3	49
31	16.5	-
-	16.7	50
32	17.0	51
-	17.3	52
33	17.5	-
-	17.7	53
34	18.0	54
-	18.3	55
35	18.5	-
-	18.7	56
36	19.0	57
-	19.3	58
37	19.5	-
-	19.7	59
38	20.0	60

※ 上記換算表は一般的なものですので、弊社製品では若干の違いがある場合がございます。